
Oxidative Behavior Materials by Thermal Analytical Techniques

Alan T. Riga and Gerald H. Patterson, editors (I) STP 1326 **STP 1326**

Oxidative Behavior of Materials by Thermal Analytical Techniques

Alan T. Riga and Gerald H. Patterson, editors

ASTM Publication Code Number (PCN): 04-013260-50



Printed in the U.S.A.

Library of Congress Cataloging-in-Publication Data

Oxidative behavior of materials by thermal analytical techniques / Alan T. Riga and Gerald H. Patterson, editors. p. cm.—(STP : 1326) "ASTM publication code number (PCN): 04-013260-50." Includes bibliographical references and index. ISBN 0-8031-2483-X 1. Thermal analysis—Congresses. 2. Oxidation—Congresses. ł. Riga, Alan T. II. Patterson, Gerald H., 1937– III. Series: ASTM special technical publication ; 1326. QD117.T4095 1997 620.1'173—DC21 97-39992 CIP

Copyright © 1997 AMERICAN SOCIETY FOR TESTING AND MATERIALS, West Conshohocken, PA. All rights reserved. This material may not be reproduced or copied, in whole or in part, in any printed, mechanical, electronic, film, or other distribution and storage media, without the written consent of the publisher.

Photocopy Rights

Authorization to photocopy items for internal, personal, or educational classroom use, or the internal, personal, or educational classroom use of specific clients, is granted by the American Society for Testing and Materials (ASTM) provided that the appropriate fee is paid to the Copyright Clearance Center, 222 Rosewood Dr., Danvers, MA 01923; Tel: (508) 750-8400, online: http://www.copyright.com/.

Peer Review Policy

Each paper published in this volume was evaluated by two peer reviewers and at least one of the editors. The authors addressed all of the reviewers' comments to the satisfaction of both the technical editor(s) and the ASTM Committee on Publications.

To make technical information available as quickly as possible, the peer-reviewed papers in this publication were prepared "camera-ready" as submitted by the authors.

The quality of the papers in this publication reflects not only the obvious efforts of the authors and the technical editor(s), but also the work of these peer reviewers. The ASTM Committee on Publications acknowledges with appreciation their dedication and contribution to time and effort on behalf of ASTM.

Printed in Philadelphia, PA October, 1997

Foreword

This publication, Oxidative Behavior of Materials by Thermal Analytical Techniques, contains papers presented at the symposium of the same name, held in New Orleans, Louisiana, on 21-22 Nov. 1996. The symposium was sponsored by ASTM Committee E37 on Thermal Methods. Alan Riga and Gerald Patterson served as chairpersons of the symposium and are editors of the resulting publication.

Contents

Overview	vii
Oxidative Behavior of Polymers and Petrochemicals	
Keynote Address	
Oxidative Induction Time—A Review of DSC Experimental Effects— R. L. BLAINE, C. J. LUNDGREN, AND M. B. HARRIS	3
Keynote Address	
Oxidation Studies of Lubricants from DORNTE to Microreactors—J. M. PEREZ	16
Thermal Studies on the Combustion Process of Treated and Untreated Cellulose—P. AGGARWAL AND D. DOLLIMORE	29
Further Studies on Oxidative Induction Test on Medical Polymers—L. woo, S. Y. DING, M. T. K. LING, AND S. P. WESTPHAL	44
The Use of Thermal Analysis to Assess Oxidative Damage in Polyolefins— A. R. HORROCKS, J. MWILA, AND M. LIU	58
Evaluation of the Oxidation Behavior of Polyethylene Geomembranes Using Oxidative Induction Time Tests—Y. G. HSUAN AND Z. GUAN	76
Thermal Oxidative Behavior of Readily Available Reference Polymers A. T. RIGA, R. M. COLLINS, AND G. H. PATTERSON	91
Thermoanalytical Approach for Differentiation Between Surface and Bulk Oxidation for Polyethylene Films —I. DOLGOPOLSKY, Y. I. GUDIMENKO, AND J. I. KLEIMAN	102
Thermo-Oxidative Degradation of Irradiated Ethylene Tetrafluoroethylens C. J. WOLF, S. C. HAGER, AND N. P. DEPKE	116
Thermal Oxidative Stability of Fiber-Reinforced Epoxy Composite Materials—M. F. FLESZAR	128
Laboratory Simulation of Oxidative Aging of Asphalt—G. M. MEMON AND B. H. CHOLLAR	138
INSTRUMENTAL TECHNIQUESOXIDATION INDUCTION PROCESSES	
Development of a Standard Test Method for Determining Oxidative	

 Induction Time of Hydrocarbons by Differential Scanning Calorimetry

 and Pressure Differential Scanning—A. T. RIGA AND G. H. PATTERSON
 11

Comparative Techniques for Oxidative Induction Time (OIT) Testing—	
R. B. CASSEL, A. W. SALAMON, G. CURRAN, AND A. T. RIGA	164
Estimation of Bias in the Oxidative Induction Time Measurement by Pressure	
DSC—S. M. MARCUS AND R. L. BLAINE	172
Oxidative Induction Time (OIT) Determinations of Polyethylenes: Influences of Temperature, Pressure, and Crucible Materials on the Result—	
R. TRUTTMANN, K. SCHIANO, AND R. RIESEN	184
A Proposed Reference Material for Oxidative Induction Time by Differential	
Scanning Calorimetry—R. L. BLAINE AND M. B. HARRIS	193
Abbreviated Oxidative Induction Time Method Without Heating/Cooling and	
Gas Switching—r. riesen, r. truttmann, and k. schiano	205
Oxidation as a Contamination in Measurements Involving Self-Generated	
AtmospheresR. J. SEYLER AND R. MOODY	212

Overview

This Special Technical Publication represents a compliation of presentations from an international symposium addressing the Oxidative Behavior of Materials by Thermal Analytical Techniques which was held 20-21 Nov. 1996 in New Orleans, Louisiana. The symposium and this subsequent publication examine the new thermal analytical techniques describing the physical properties and oxidative degradative behavior of polymers, lubricants, and petrochemicals. Historical reviews, oxidation mechanisms, new test methods, unique techniques, robotic methods, new reference standards, and bias considerations form the basis of this publication.

It was generally agreed, but certainly highlighted by Roger Blaine, TAI, keynote speaker, that the factors affecting the oxidation induction time (OIT) were: the isothermal temperature; pan type, metallurgy, and shape; and pressure and oxygen flow rate. Polymer oxidation was emphasized in a number of papers. Professor Joseph Perez, Penn State University, keynote speaker, focused on the measurement of oxidation in formulated passenger car and diesel engine oils. Professor Perez discussed the evolution of oxidation systems from the complex Dornte type systems of the 1940s to the current use of microreactors, for example, DSC, TGA, and Klaus Penn State Microreactor (PSMO). He stressed the variables such as: metal surfaces that affect the rate of oxidation or thermal decomposition, temperature effects on the related thermal and oxidative processes, effects of volatility on vapor phase versus liquid phase oxidation, oxygen diffusion rate limitations, and additive effectiveness. Some of the differences in the thermal and oxidative behavior in bulk systems (engine tests) and microsystems (DSC, TGA, PSMO) were discussed.

Alan Riga, Lubrizol, reviewed the recently approved and soon to be published Standard Test Method for Determing OIT of Hydrocarbons by DSC/PDSC. An interlaboratory study using this DSC/PDSC method was reported by Blaine and Riga, with ASTM Reference C, a diluted fully formulated engine oil from Lubrizol and ASTM Reference D, a polyethylene film from TAI.

A successful robotic DSC evaluation by M. Kelsey, Mettler-Toledo, clearly differentiated References C and D at 195°C and one atmosphere of oxygen. This abbreviated OIT method without heating/cooling (from room temperature) or gas switching reduces the experimental time, avoids a nitrogen purge, and gives a slightly better reproducibility. The abbreviated OIT method suggests a possible revision to the DSC ASTM E37 test protocol. A question raised by a number of participants was "Is it possible that some PDSC oxidative testing can be replaced with robotic DSC?"

A modified Arrehenius model was used to predict OIT of polymers. An adaptation of this model for evaluation of the oxidative stability of oils was discussed.

A new scanning DTA/TGA oxidation test was presented. This DTA method is based on air oxidation and defines an oxidation temperature. A good correlation was observed between the isothermal PDSC OIT in oxygen of eight readily available olefin reference polymers and the oxidation temperature. The DTA air oxidation method as well as the PDSC oxygen method use the melting temperatures and heat of fusion of polyethylenes and polypropylenes to verify the temperature and heat calibrations.

A unigue approach to the oxidation process was presented by Rick Seyler (Kodak). This paper considered oxidation as a negative factor in the determination of vapor pressure by

DSC in nitrogen. The observed exotherms in DSC measurements were associated with partial oxidation of the chemical specimen from residual air in the DSC specimen container. An application of this method is to study vapor phase oxidation of organics.

Other applications of DSC or PDSC oxidation tests included: radiation-damaged polyethylene, cellulose, thin film oxidation, medical poymers, and asphalts.

The attendees agreed that this symposium was rewarding and much knowledge was gained. They commented that the prsentations focused on areas that have not been previously discussed. Symposium cochairman, Alan Riga, suggested that another meeting be organized in two years at a future NATAS conference. The speakers and attendees agreed.

Rick Seyler, Kodak, summarized the symposium in the form of the following poem:

OIT BLUES

A SYMPOSIUM REVIEW

By R. J. Seyler

Don't you know its not thermodynamic? Rather, OIT is really quite kinetic!

So when you do your OIT It most likely won't be like me!

With so many ways to do the test Precision within-lab will be the best.

When finally on a protocol we agreed Eaxct control of temperature we'll need.

So now proceed as best we can Just make sure its an aluminum pan.

Beware when other metal is present Copper catalysis is a serious contaminant.

Other conditions like area, weight, and flow Their influence we will need to know.

To switch the purge gas it was shown One need not, for the BIAS is now known.

A last condition we yet must know What point in time do we declare zero?

When all our tests we run one way Which value of OIT do we report today?

Extrapolated onset, threshold, or the peak 'Tis a single number of value that we seek.

We need to realize that induction time May have different significance for yours and mine.

Unless our sample is without formulation OIT will address the type of stabilization!

We have tried our best to listen well To learn all the messages you had to tell. But when all is said and done We should yet be troubled by just this one.

That from accelerated aging over which we toiled The egg we hatched may have been Hard Boiled!

Alan T. Riga The Lubrizol Corporation Wickliffe, OH; Symposium cochairman and coeditor

Gerald H. Patterson The Lubrizol Corporation Wickliffe, OH; Symposium cochairman and coeditor

ISBN 0-8031-2483-×