

SUBJECT INDEX

A

Adsorption method: 30, 37
 Alcohol sulfates: 4, 6, 7, 10, 11,
 17, 25, 43, 51, 69, 85
 Alkyl aryl sulfonate: 27, 71, 72,
 73, 79, 85, 94
 Alkyl benzene sulfonate: 27
 Alkyl sulfate: 4

B

Brilliant Aviroil: 10

C

Cationic agents: 18, 23, 26, 32,
 40, 42, 44, 45, 48, 52, 54, 58,
 68, 70, 74, 87, 90, 93
 Colorimetric methods: 26, 27, 34,
 42, 44, 65, 69, 76, 87, 93
 Conductimetric method: 24

D

Drop number: 39

F

Functional groups, aid in
 analysis: 78, 82

G

Gardinol: 3, 8, 10

I

Igepon: 3, 8, 10
 Infrared methods: 82

M

Moisture: 81, 94

N

Nekal: 8

Nonionic detergents: 20, 60, 74,
 79, 86, 92

P

Polarographic adsorption analysis:
 28

Q

Qualitative methods: 7, 8, 12, 19,
 36, 38, 46, 47, 56, 57, 64, 66,
 67, 77, 79, 80, 83, 86, 89, 92
 Quantitative methods: 1-4, 6, 9-11
 13-18 21-27, 29-35, 37, 39, 40-
 45, 48-52, 54, 55, 58-63, 65,
 68-76, 78, 80, 81, 84, 87, 88,
 90, 93-96
 Quaternary ammonium compounds: 26,
 32, 40, 42, 44, 45, 48, 52, 54,
 58, 68, 70, 74, 87, 90, 93

S

Soap-synthetic detergent mixtures:
 29
 Sodium sulfate: 24
 Sulfated alcohols: 4, 6, 7, 10, 11,
 17, 25, 43, 51, 69, 85
 Sulfonated oils: 2, 3, 8, 10, 14,
 15, 61
 Sulfonic acid salts: 9, 27, 33,
 34
 Sulfur: 13
 Sulfuric anhydride: 1, 2, 3, 4,
 15, 16
 Surface tension method: 31, 35

T

Titration with cationic agent
 (or vice versa): 35, 48, 50,
 51, 62, 63, 70, 75, 84, 95,
 96
 Turbidity methods: 48, 52
 Turkey-red oils: 2, 3, 8, 10, 14,
 15, 61