

## Subject Index

### A

- Abrasion, liner, 177
- Acetabular components, 177
  - cups, 60, 199
  - cups, cementless, 21
- ASTM standards, 177

### B

- Ball/stem interface, 127
- Bone support test, partial, 137
- Bore and cone interface, 189
- Burnishing, 177
- Burst strength, 127

### C

- Cantilever test, 137
- Ceramic femoral head, 69
- Chromium, 45
- Chromium orthophosphate, 33
- Cobalt, 45, 104
- Cobalt-chrome alloy, 33, 60, 104, 127, 146, 157
- Cobalt chromium molybdenum alloy, 45, 146, 189, 211
- Compressive stress, 94
- Corrosion, 4, 45
  - bore and cone interface, 189
  - crevice, 211
  - electrochemical changes, 157
  - femoral head and neck taper interface, 104
  - fretting, 146, 211
- Morse taper, 114
- product characterization, 33

### D

- Deformation, 85, 127
- Diameter/shell inside diameter, 60
- Disassembly force, Morse taper, 114

- Displacement vs. load information, 94
- Distraction resistance, 199

### E

- Energy dispersive X-ray analysis, 33

### F

- Fatigue characteristics, 226
- Fatigue integrity test, 137
- Fatigue loading, 104, 199
- Fatigue stem testing, 226
- Fatigue testing, 211
- Fatigue, torsional, 177
- Femoral component, 114, 137
- Femoral heads, 69
  - modular, 33
- Femoral stem fatigue, 226
- Fibroblast, 21
- Fourier transform infrared spectroscopy, 33
- Fretting, 33, 157
  - bore and cone interface, 189
  - corrosion fatigue, 146
  - corrosion, implant interface, 211
- intercomponent, 226
- modular connection, 137
- taper, crevice, 45

### G

- Glenoid prosthesis, 5

### G

- Head-neck taper interface, 104
- Head-stem combination, 5
- Head-taper joints, 146
- Hip joint head, 127
- Hip joint replacement, 127

Hip replacement, 33, 45, 211,  
    226  
    acetabular cups, 21, 60  
    ceramic femoral head  
        marking, 69  
    corrosion testing, 157  
    fatigue loading, 199  
    fretting corrosion, 189  
    interface fretting assessment,  
        104  
    load bearing capacity, 127  
    modular joint corrosion, 146  
    Morse taper function, 114  
    torsional fatigue, acetabular  
        components, 177  
Humeral prosthesis, 5  
Hybridization, 21

**I**

Immunohistochemistry, 21  
Impaction force, Morse taper,  
    114  
International Organization for  
    Standardization (ISO), 177

**K**

Knee systems, 85, 94, 137

**L**

Laser beam engraving, 69  
Liner lock mechanism, 60, 177  
Liner micromotion, 177  
Liner/shell distraction resistance,  
    199  
Load bearing capability, 127  
Loading conditions, 226  
Loading, 226  
    Loading, fatigue, 85, 94  
    Load magnitude effects, 157  
    Load resistance, static, 127  
Locking mechanism, 94  
Locking strength, 85, 94  
Loosening, aseptic, 21, 60

**M**

Machining, mechanical, 69  
Marking, ceramic femoral head,  
    69

Modular junctions, 33  
Morse taper, 104  
    assessment, 114

**N**

Neck extension, 189  
Neck taper, 104

**O**

Open circuit potential, 157  
Osteolysis, 21, 33, 60  
Oxide film, 157  
    fracture, 45

**P**

Particulate debris, 33  
Physiologic loading, 85  
Polyethylene, 60, 177  
    insert failure, 5  
    liner, 199  
    wear, 21, 33

**R**

Ringer's solution, 104, 146

**S**

Shear load, 85  
Shear stress, 94  
Shoulder arthroplasty, 5  
Stability, modular connection,  
    94  
Stem extension test, 137  
Stress concentrations, 69  
Structural fatigue strength, 226  
Surface abrasion, 60  
Surface morphology, 146

**T**

Taper assessment, 104, 114,  
    146, 157  
    corrosion, 5  
    distraction strength, 114  
    fit, 211  
    fretting crevice corrosion, 45  
    interface, 104  
    surface structure, metal, 127

Tibial components, 85, 94  
Tibial stem extension, 137  
Titanium alloys, 33, 104, 157  
  heads, 45  
  stems, 45, 189, 211  
  tapers, 127, 146  
Torsional fatigue, 177

**W**  
Wear, 33, 177

Wear debris, modular hip implants, 114

**X**

X-ray diffraction analysis, 33

**Z**

Zirconia ceramic femoral head,  
  146