

Index

A

- Albumin, 105
- Anderson, J. M., 221, 251, 321
- Aneurysm clips, 147
- Ankle joints, 424
- Anodic polarization, 15
- Arterial prostheses, 294
- Ashford, R., 373
- Asphahani, A. I., 30
- ASTM
 - Committee F-4 on Medical and Surgical Materials and Devices, 429
 - Practice for Assessment of Compatibility of Metallic Materials for Surgical Implants with Respect to Effect of Materials on Tissue (F 361-80), 158
 - Recommended Practice for Constant Amplitude Axial Fatigue Tests of Metallic Materials (E 466-76), 32
 - Recommended Practice for Experimental Testing for Biological Compatibility of Metals for Surgical Implants (F 361-72), 226
 - Recommended Practice for Performing Stress-Corrosion Cracking Tests in a Boiling Magnesium Chloride Solution (G 36-73) (1979), 77
 - Recommended Practice for Standard Reference Method for Making Potentiostatic and Potentiodynamic Anodic Polarization Measurements (G 5-78), 33, 79, 121
 - Recommended Practice for Static Bend Testing of Nail Plates (F 384-73), 93
 - Recommended Practice for Surface Preparation and Marking of Metallic Implants (F 86-76), 139, 353
 - Specification for Cast Cobalt-Chromium-Molybdenum Alloy for Surgical Implant Applications (F 75-76), 32, 410
 - Specification for Hip Nail—Jewett Type (F 369-73), 93
 - Specification for Stainless Steel Bar and Wire for Surgical Implants (F 55-76), 59, 74
 - Specification for Stainless Steel Bars and Wire for Surgical Implants (Special Quality) (F 138-76), 59, 108, 197, 210
 - Specification for Stainless Steel Sheet and Strip for Surgical Implants (F 56-78), 119
 - Specification for Titanium 6Al-4V ELI Alloy for Surgical Implant Applications (F 136-79), 74, 93, 108
 - Specification for Unalloyed Titanium for Surgical Implant Applications (F 67-77), 108

- Specification for Wrought Cobalt-Nickel-Chromium-Molybdenum Alloy for Surgical Implant Applications (F 562-78), 108, 197
- Tension Testing of Metallic Materials (E 8-79), 32
- Test for Nickel in Water (D 1886-77), 107
- Test for Rubber Deterioration—Crack Growth (D 813-59), 280
- Test for Strength Properties of Adhesives in Shear by Tension Loading in the Temperature Range from -267.8 to -55°C (-450 to -67°F) (D 2557-72), 184
- Test for Tensile Properties of Plastics by Use of Microtensile Specimens (D 1708-66), 324
- Test Method for Plane-Strain Fracture Toughness of Metallic Materials (E 399-78a), 76
- Test Method for Rubber Property—Tear Resistance (D 624-73), 280
- Test Methods for Flexural Properties of Unreinforced and Reinforced Plastics and Electrical Insulating Materials (D 790-81), 383
- Test Methods for Rubber Properties in Tension (D 412-75), 283
- Tests for Flexural Properties of Plastics and Electrical Insulating Materials (D 790-71), 344
- Atomic absorption spectrophotometry, 223
- B**
- Ball-joint rods, 160
- Bartel, D. L., 415
- Bearden, L. J., 208
- Bend tests, 61
- Biocompatibility of materials, 30, 208, 223, 251, 278
- Biodegradable materials, 340
- Black, J., 28, 71, 223, 265
- Blais, P., 294
- Bone analysis, 387
- Brown, C. A., 160
- Brown, S. A., 105, 195
- Bundy, K. J., 73
- Burstein, A. H., 415
- C**
- Cage implant system, 251
- Carbon, 274
- Carbon fibers, 340
- Carboxylic acids, 294
- Casper, R. A., 340
- Cells, 208
- Cement stresses, 373
- Ceramic fibers, 340
- Ceramics, 369
- Charnley hip joint, 417
- Chromium alloys
- Co-Cr-Mo alloys, 11, 251, 373
 - Co-Ni-Cr-Mo alloys, 105
- Corrosion reaction, 5
- Grommets, 274
- Metabolism, 242
- Properties, 30
- Systemic effects, 223
- Tissue sensitivity, 209
- Zimaloy, 11
- Clark, G. C. F., 106
- Clarke, I. C., 351, 373
- Coatings, 6, 31, 179
- Cobalt alloys
- Aneurysm clips, 147
 - Cast, 12, 30
 - Co-Cr-Mo alloys, 351, 373
 - Co-Ni-Cr-Mo alloys, 418
 - Corrosion, 427
 - Grommets, 274
 - Properties, 30

Tissue sensitivity, 209
 Zimaloy, 11
 Composite materials, 340
 Consensus development and standards,
 434

Corrosion

Forms of, 6
 Crevice attack, 6
 Electrochemical, 168
 Fatigue, 6, 30, 93, 179
 Fretting, 93, 105
 Oxidation effects, 28
 Pitting, 6, 117, 147, 427
 Stress-corrosion cracking, 6, 57,
 73, 147
 Potential shifts, 11
 Products, 195, 223
 Reactions, 5
 Coulometers, 168
 Cowsar, D. R., 340
 Cracking
 Fatigue, 93
 Polyurethanes, 308, 322
 Silicone, 278
 Stress corrosion, 12, 30, 57, 73, 147
 Crevice potential, 11

D

deGroot Swanson, G., 267
 Dental materials, 168
 Depolymerization, 294
 Desai, V. H., 73
 Design standards, 429
 Dimensional stability, 294
 Donald, G. D., 160
 Dujovny, M., 147
 Dunn, R. L., 340

E

Ebramzadeh, E., 373
 Edwards, B. J., 11
 Elastomers, silicone, 267, 279, 424
 Elbow joints, 424

Electrochemical corrosion, 168
 Electropolishing, 136
 External fixation, 160

F

Fatigue of materials, 5, 415, 435
 Aneurysm clips, 147
 Ball-joint rods, 160
 Cement, 373
 Co-Cr-Mo alloys, 11, 30
 Co-Ni-Cr-Mo alloys, 418
 Fretting, 93
 Polyethylene, 351
 Polyester, 294
 Silicone, 267
 316L stainless steel, 57, 73, 93, 136,
 418
 Ti-6Al-4V ELI, 73, 93, 179
 Tricalcium phosphate ceramic, 369
 Femoral prostheses, 373
 Femoral stem fracture, 403
 Fiber composites, 340
 Fibrous capsule, 267
 Film, passive, 73
 Finger joints, 278, 424
 Fissuring, 308, 322
 Flexible hinges, 278
 Flexural durability, 278
 Flow propagation resistance, 278
 Fracture mechanics (*see also* Crack-
 ing), 73, 415
 Fraker, A. C., 106, 179, 432
 Fretting
 Co-Cr-Mo alloys, 31
 MP35N, 105
 316L stainless steel, 93, 105
 Ti-6Al-4V alloy, 93, 105
 Titanium, 105
 Frisch, E. E., 278

G

Galvanic activity, 147
 Garton, A., 294

Georgette, R., 321
 Giant-cell response, 267, 425
 Gilbert, J. L., 93
 Gold alloys, 168
 Greene, N. D., 5, 55, 71, 88, 115
 Grommets, 274
 Guidoin, R., 294
 Gunasekera, K. R., 294
 Gustavson, L., 56, 190

H

Hahn, H., 179
 Hank's solution, 117
 Hickl, J., 30
 Hiltner, A., 251
 Hip prostheses, 30, 93, 351
 Hood, A. C., 56
 Hood, R. W., 403
 Horowitz, E., 265, 320, 427
 Horvath, N., 322
 Host tolerance, 267
 Hydrogen embrittlement, 11
 Hydrogen uptake, 13

I

Implant retrieval, 403, 415
 Incentives, marketplace, 434
 Infection, 195
 Inflammatory response, 251
 Infrared spectroscopy, 294
 Iron, 208, 223, 237
 Irving, C. C., 136
 Irwin, G. R., 74

J

Joint components, 403, 415
 Joint replacement, 373, 415
 Joint simulator wear tests, 351

K

Kafesjian, R., 207
 Kies, J. A., 74

King, M. W., 291, 294
 Knee joints, 422
 Kossovsky, N., 147
 Kossowsky, R., 147
 Kumar, P., 30

L

Langley, N. R., 278
 Lare, P. J., 179
 Lautenschluger, E. P., 223
 Lawley, A., 30
 Lemons, J., 25, 89, 178, 208, 369
 Liability, product, 434
 Lipid absorption, 278
 Loosening of hip prostheses, 373
 Louthan, M. R., Jr., 11
 Lucas, L. C., 208

M

McArthur, W. A., 308
 McCoy, K. D., 322
 McKellop, H. A., 351
 Malpractice, medical, 434
 Mandible repair, 340
 Manganese, 208
 Marchant, R. E., 251
 Marek, M., 168
 Maupin, B. K., 267
 Maxillofacial surgery, 340
 Mayor, M. B., 429
 Merritt, K., 105, 195
 Metal alloys (*see also* Steel, Titanium)
 Chromium-nickel, 5
 Corrosion products, 195
 Corrosion reaction, 5
 Tissue sensitivity, 195
 Transport in the body, 203
 Zimaloy (Co-Cr-Mo), 11
 Metal-cell complexes, 195
 Metal-protein complexes, 195
 Metal salts, 195, 208
 Mina-Araghi, M., 373

Molybdenum alloys
 Co-Cr-Mo alloys, 11, 351
 Co-Ni-Cr-Mo alloys, 105
 Properties, 30
 Tissue response, 209
 Zimaloy, 11
 Morin, C. R., 57
 Mueller hip joint, 417, 419

N

N-formyl-Met-Leu-Phe (MLP), 251
 Nails, 93
 Nalbandian, R. M., 267
 Nickel alloys
 Co-Ni-Cr-Mo alloys, 105
 Corrosion reaction, 5
 Tissue response, 209

O

Ogundele, G. I., 117, 178, 206, 321
 Olson, R. W., 322
 Ordway, F., 179
 Organ accumulations, 223
 Oxidation (*see also* Corrosion), 308
 Oxide layer, 136

P

Pacing leads, 308, 322
 Packer, K. F., 57
 Parins, D., 264, 323
 Performance standards, 429, 434
 Piehler, H. R., 56, 71, 93, 292, 434
 Pitting, 7
 Co-Cr-Mo alloys, 31, 39
 Steel alloys, 117, 427
 Plasma spraying, 179
 Plastics (*see* Polymeric materials)
 Plates, 93, 105
 Polarization, 8, 11, 73, 117
 Polishing, 136
 Polymeric materials
 Poly(ethylene terephthalate), 294

Poly (2-hydroxyethyl-L-glutamine)
 (PHEG), 251
 Polyethylene, 274, 351, 403
 Polyesters, 294, 340
 Polyurethanes, 308, 322
 Resins, 340
 Silicone, 266, 278
 Portnoff, M. A., 93
 Powder metallurgy, 11, 30
 Product liability, 434
 Proplast, 274
 Prosthetic hips, 351
 Protein albumin, 105

R

Radiolucent zones, 373
 Regulation
 Medical device, 434
 Social, 434
 Reinforced plastics, 340
 Resins (*see also* Polymeric materials),
 185, 340
 Retrieval analysis, 403, 415
 Risk assessment and management, 434
 Rods, 160
 Rowe, R. H., Jr., 179

S

Salts, tissue response to, 195, 208
 Scaffolds, porous, 369
 Screws, 105, 108, 195
 Seligson, D., 160
 Semlitsch, M., 56, 246, 277, 292, 414,
 427
 Sensitivity of tissues, 195
 Serum, 105
 Sheehan, J. P., 57
 Shetty, H. R., 71, 115
 Silicone, 267, 278, 424
 Silver alloys, 168
 Sintering, 183
 Sisson, R. D., Jr., 11

- Slotter, L. E., 93
 Smith, G. K., 223
 Smith, T., 339
 Speck, K. M., 106
 Spires, W. P., Jr., 190
 Standards, performance, 429, 434
 Standards writing, 429
 Starkebaum, W., 414
 Steel
 Stainless
 Corrosion products, 195
 Effect of crevices on potential shifts, 13
 Effects on cells, 195, 208
 Grommets, 274
 Hip prosthesis loosening, 373
 Pitting, 427
 17-7PH stainless, 147
 18Cr-8Ni, 6
 18Cr-8Ni-2Mo, 6
 301 stainless, 147
 303 stainless, 147
 304 stainless, 58, 85-86
 316L stainless
 Aneurysm clips, 147
 Corrosion, 37, 117, 223
 Effects on cells, 208
 Electropolishing, 136
 Fatigue, 418
 Fretting corrosion, 105
 Stress corrosion cracking, 57, 73
 1018 carbon, 8, 9
 Strain rate, 57
 Stress
 Cracking, 57, 308, 322
 In cement, 373
 Residual, 320
 Stress corrosion cracking, 7, 11, 30, 73
 Sugie, T., 251
 Sulfur, 168
 Surface fissuring, 308, 322
 Surgical devices, 429
 Swanson, A. B., 267
 Szycher, M., 308

 T
 Tafel extrapolation, 8
 Tarnishing, 168
 Tear propagation strength, 278
 Temperature effects, 128, 132
 Tension tests
 Co-Cr-Mo alloys, 11, 30
 316L stainless steel, 60
 Titanium alloy, 183
 Tissue response, 195, 208
 Titanium
 Corrosion, 13, 105, 427
 Grommets, 274
 Sponge, 179
 Ti-6Al-4V alloy, 93, 105, 179, 351, 373
 Toe joints, 278
 Total-hip replacement, 373
 Transmission electron microscopy, 208
 Trapezoidal-28 hip joint, 417
 Tricalcium phosphate ceramic, 369

 U
 Ultrastructures of cells, 208

 V
 VanOrden, A. C., 264, 427
 Vascular diseases, 294
 Vascular surgery, 147
 Vegdahl, E. J., 93

 W
 Wear
 Clinical, 351
 Debris, 267, 403
 Metal, 422
 Polyethylene, 351, 403, 415
 Weber, M. J., 93

White, W. E., 117
 Williams, D. F., 106
 Wound healing, 340
 Wright, T. M., 415

X

X-ray diffraction, 147

Y

Yapp, R., 266

Z

Zimaloy, 11